

# RetCam™ Review Software with Patient Database Application



Enables remote evaluation of digital images and provides patient database application for image storage and retrieval.

#### *Enhanced Telemedicine Utility:*

Now RetCam Digital Images can be transferred directly to your networked server or physician specialist for remote review and diagnosis. The efficient transfer of images enables timely evaluation of critical ophthalmic disorders and allows in-depth analysis and comparison of patient images with the RetCam Review Software.

#### *The RetCam Review Software with Patient Database Application:*

- Enables the transport of images to any networked server, pc, or RetCam System
- Allows remote evaluation of images at any time by one or more designated reviewers
- Includes advanced image processing and image compare operations
- Enhances data back up and security of patient records



#### *Compare Screen:*

- Side-by-side longitudinal image comparison
- Pre and post treatment
- Compare with other sessions or other patients



#### *Image Processing:*

- Adjustments for optimal image analysis
- Apply viewing enhancements to single image or complete clinical session
- Export images with optional patient demographics



**Clarity**  
medical systems



## Technical Specifications

The RetCam Review Software/Patient Database Application Package can be installed in any computer that fits the following specifications:

Operating System:	Windows® XP Professional Service Pack 2 or higher
Processor Speed:	Minimum Pentium® M 2.0 GHz
Hard Disk Storage Space:	Minimum 80 GB
System Memory:	Minimum 1.0 GB
Video Memory:	Minimum 128 MB
Network Interface:	100 Megabit/s Ethernet transfer rate.
Optical Media Drive:	Writable/recordable optical media device



**Clarity**  
medical systems